OMRON®

# INSTRUCTION MANUAL

Pedometer with Download Capability

Model **HJ-720ITCAN**

FRANÇAIS

ENGLISH

# CONTENTS

## Before Using the Unit

Introduction .....	3
Important Safety Information.....	4
Know Your Unit .....	6
Main Unit .....	6
Components .....	7
Aerobic Step Function .....	8
Calorie/Amount of Fat Burned Display .....	8
Stride Length Measurement .....	9

## Operating Instructions

Battery Installation .....	10
Battery Replacement .....	11
Setting the Time .....	12
Setting the Weight .....	14
Setting the Stride Length.....	15
Attaching the Pedometer.....	16
Changing the Display .....	18
Using the Event Function.....	18
Using the Memory Function .....	19
Using the System Reset Function .....	21
Downloading Data to a PC.....	22

## Care and Maintenance

Care and Maintenance.....	23
Troubleshooting Tips .....	24
Specifications .....	25
FCC Statement .....	26
Warranty Information.....	27

# INTRODUCTION

Thank you for purchasing the Omron® HJ-720ITCAN Pedometer with download capability.

**Fill-in for future reference**

**DATE PURCHASED:** \_\_\_\_\_

**Staple your purchase receipt here.**

A pedometer is a great motivational tool to track your distance, the number of steps you take and how many calories you burn. The convenient memory function supports the management of daily walking. You can easily set and reach your daily goals. This pedometer allows the user to place the pedometer in a pocket, bag or clip it to a belt when walking for added flexibility to count steps.

Your HJ-720ITCAN comes with the following components:

- Pedometer
- Screwdriver
- Strap
- Clip
- Holder
- Battery (CR2032)
- Instruction Manual

## SAVE THESE INSTRUCTIONS

# IMPORTANT SAFETY INFORMATION

*To assure the correct use of the product basic, safety measures should always be followed including the warnings and cautions listed in this instruction manual.*

## SAFETY SYMBOLS USED IN THIS INSTRUCTION MANUAL

 <b>WARNING</b>	Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.
 <b>CAUTION</b>	Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property.

## OPERATING THE DEVICE

- ⚠ Keep the unit out of reach of young children. If a young child swallows the battery, battery cover, or screw immediately consult with a doctor.
  
- ⚠ Read all information in the instruction book and any other literature included in the box before using the unit.
  
- ⚠ Operate the unit only as intended. Do not use for any other purpose.
  
- ⚠ Properly dispose of used battery. Do not throw the battery into fire. The battery may explode.
  
- ⚠ Do not swing the pedometer by holding the strap. This may cause injury.


# IMPORTANT SAFETY INFORMATION

## CARE AND MAINTENANCE


- ⚠ Do not disassemble or modify the unit. Changes or modifications not approved by Omron Healthcare will void the user warranty.
- ⚠ Consult your physician or healthcare provider before beginning an exercise program.

# KNOW YOUR UNIT

## Main Unit

## Back of the main unit


# KNOW YOUR UNIT

## Components


### Holder

Use this holder when attaching the unit to your belt or the top of your pants.


- 1) Attach the main unit to the holder.


- 2) Fasten the holder clip to your pants.

### How to attach the strap to the main unit


### How to connect the strap and the clip

### How to release from a clip and fasten with a clip

To fasten

To release

# AEROBIC STEP FUNCTION

The unit displays the total number of aerobic steps in a day that satisfies the following two conditions: (The minutes you walked as the number of aerobic steps is also displayed.)

- **Walk more than 60 steps per minute**
- **Walk for more than 10 minutes continuously**

Example: If you walk for 20 minutes at the speed of 120 steps per minute, the number of steady aerobic is 2400.

**Note:** Taking a rest for less than 1 minute during continuous walking after 10 minutes is considered as continuous walking.

# CALORIE/AMOUNT OF FAT BURNED DISPLAY

The calorie and amount of fat burned display indicates a calorie consumed by walking and length from the weight, stride distance, number of steps, and walking speed.

The calorie and amount of fat burned display increment is one minute.

If you want to know the target calorie value, multiply the displayed calorie by the coefficient in the table shown right.

## Example of correction:

For a woman in her 50's with the calorie display of 300, the consumed calorie is calculated as follows.

$$300 \times 0.95 = 285 \text{ kcal}$$

## Correction coefficient for calorie display

	Man	Woman
20's	1.10	1.07
30's	1.05	1.01
40's	1.02	.97
50's	1.01	.95
60's	1.00	.95
70's	.98	.96
80's	.95	.95


# STRIDE LENGTH MEASUREMENT

Measure your stride length before beginning to set the pedometer.

## How to measure your stride length

1. Walk 10 steps with your normal stride as illustrated below.

2. Measure the distance from START to END in centimeters.

3. Calculate your stride length by dividing the total distance by 10.

### Example:


Total distance = 800 cm

$800 \text{ cm} \div 10 = 80 \text{ cm}$


# BATTERY INSTALLATION

## 1. Remove the screw on the battery cover located on the back of the main unit.

Use the supplied screwdriver or use a commercially available small screwdriver to remove the screw on the battery cover.


## 2. Open the battery cover.


## 3. Insert the battery.

Place the positive (+) side of the battery (CR2032) upward. Insert the battery in the direction of the arrow as illustrated.

## 4. Close the battery cover.

Insert the projected portion of the battery cover first.


## BATTERY INSTALLATION

### 5. Tighten the screw on the battery cover.

#### NOTES:

- The time setting display will appear on the screen for approximately 5 minutes after the battery is inserted.
- After installing the battery, set time of day, weight and stride length. (Refer to “Setting the Time” on page 12)

## BATTERY REPLACEMENT

When the battery is removed from the unit, the time of day and the data stored in the memory for the current hour will be deleted. The weight, stride length setting and the data stored in the memory up to the previous hour will not be deleted.

The current day is automatically stored in the memory when the time reaches 0:00. Do not replace the battery at 0:00. For example: If you remove the battery at 23:59 and replace the battery at 0:01 the data for the previous day will be stored as the current day.

Low Battery Indicator


When the Low Battery Indicator appears on the display, remove the old battery.

#### NOTE:

When the battery indicator starts to blink, the unit stops counting and the time is displayed as “-:-:-”.

Take the battery out with a thin stick that is not easily breakable.

- DO NOT use a pair of metal tweezers or a screwdriver.
- Replace with a new CR2032 battery.
- Dispose of a battery according to applicable local regulations.

# SETTING THE TIME

**Setting range: 0:00 to 23:59**

## 1. Install the battery.

(Refer to "Battery Installation" on page 10)

The hour flashes on the display.


If the MEMO/▲ or SET Button is not pressed for 5 minutes, complete setting and display the number of steps and time.

## 2. Press the MEMO/▲ Button to advance by increments of one hour.


### NOTE:

Press the MODE/EVENT button to return the hour to 0.


## 3. Press the SET Button to set the hour.

The minutes flash on the display.

## 4. Press the MEMO/▲ Button to advance by increments of one minute.

Press and hold the button to advance the time in 10 minute increments.


### NOTE:

Press the MODE/EVENT Button to return the minute to 00.

## SETTING THE TIME

**5. Press the SET button to set the minute.**


The weight flashes on the display.

# SETTING THE WEIGHT


**Setting range: 30 to 136 kg.**

The weight flashes on the display.

- 1. Press the MEMO/▲ Button to advance by increments of 1 kg.**

Press and hold the button to advance weight in 10 kg increments.


**NOTE:**

Press the MODE/EVENT button to return the weight to 30.

- 2. Press the SET Button to set the weight.**


The stride length flashes on the display.

# SETTING THE STRIDE LENGTH


## Setting range: 30 to 120 cm

The stride length flashes on the display.

### 1. Press the MEMO/▲ Button to advance by increments of 1 cm.

Press and hold the button to advance the stride in 10 cm increments.

#### NOTE:

Press the MODE/EVENT Button to return the stride length to 40 cm.

### 2. Press the SET Button to the stride length.

The number of steps and time will display.

**All the settings are completed.**

**To change the time, weight, or stride length after the setting is completed:**

1. Press and hold the SET Button for 2 seconds.
2. Change the setting.

#### NOTE:

If setting does not need to be changed, press the set button to advance to the next setting.

# ATTACHING THE PEDOMETER

Place the unit in the upper front pocket of your jacket, in a bag you hold tightly, or attach unit to your pants or belt.

## Place the unit in your pocket.

- Do not place the unit in the back pocket of your pants.
- Attach the clip on the edge of your pocket to prevent the unit from dropping.


## Place the unit in your bag.

- Place the unit in your pocket or the partition of your bag.

## Attach the unit to your belt or the top of your pants.

- Use the supplied holder.
- Do not force the holder onto thick or wide.

### NOTE:

The display will not change until you have walked for more than 4 or more seconds.


# ATTACHING THE PEDOMETER

**NOTE: The unit may not be able to count the steps correctly in the following cases.**

**When the memory is being displayed or the unit is being set**


**When the unit is not vertical to the ground**

- When the front of the main unit is placed in the angle of less than  $60^\circ$  (shown to the right) or horizontal to the ground.

**NOTE:**

The unit can count steps even if it is slanted or placed upside-down, or if the main unit is perpendicular to the ground.

**When the main unit moves irregularly**

- When the main unit is placed in the bag that moves irregularly because it hits your foot, your belt or the top of your pants.
- When the unit hangs from your belt, the top of your pants or a bag.

**When you walk at an inconsistent pace**

- When you shuffle or wear sandals
- When you fail to walk consistently in a crowded place

**When you use the unit in the place where lots of up and down movements or vibrations take place**

- Standing up and/or sitting down movement
- Playing sports other than walking
- Ascending or descending movement at stairs or in a steep slope
- Vertical or horizontal vibration in the moving vehicle such as on the bicycle, in the automobile, trains or a bus

**When you jog or walk extremely slowly**



**NOTE:**

The unit will count steps during jogging if it is attached to your belt or the top of your pants with a holder.

## CHANGING THE DISPLAY


**Press the MODE/EVENT Button to change the display**

Each time you press the MODE/EVENT Button, the display will change.

## USING THE EVENT FUNCTION

**Press the MODE/EVENT Button for 2 seconds while current day's data is displayed.**

The times of event in current day is displayed for 2 seconds.

### NOTES:

- Follow the direction by instructor or service provider to use this function.
- Event can be recorded once per hour.
- The unit does not display previous event data.

# USING THE MEMORY FUNCTION

This unit can store up to previous 41 days and display the most recent 7 days of data.

The current day's data is automatically stored in memory when the time reaches 0:00. The display will return to 0.

## NOTE:

The time must be correctly set in the unit.

## How to view data stored in the memory


**Press the MODE/EVENT Button** until the data you want to see is displayed:

Steps

Aerobic Steps with minutes

Calories with amount of fat burned

Distance


**Press the MEMO/▲ Button.** The previous day data is displayed.

Continue to press the MEMO/▲ Button until all seven days of data is displayed.

## NOTE:

If no button is pressed for more than 1 minute, the display will return to the current day's data.

**Press the MODE/EVENT Button** to return to the current day's data.

# USING THE MEMORY FUNCTION

## Full Memory Indicator

### NOTES:

- Full memory indicator starts display when there are 35 days including current day of stored data which has not been downloaded.
- Follow the direction by instructor or service provider to download the data. (Refer to “Downloading Data to a PC” on page 22)

# USING THE SYSTEM RESET FUNCTION

Use the **SYSTEM RESET** function in the following cases:

- When the display becomes unrecognizable


**Press the SYSTEM RESET Button on the back of the main unit with a thin stick that is not easily breakable.**

The time is reset.

**NOTE:**

The data stored in the memory is not deleted.

After all characters on the displays are turned on, the initial time setting displays appears. (Refer to “Setting the Time” on page 12)


# DOWNLOADING DATA TO A PC

Use only the authorized Omron USB cable sold separately, (Model Number HHX-CABLE-USB2) with this pedometer. Use of any other cable may cause damage to the pedometer and will void the user warranty.

## 1. Open the USB Connector cover.


## 2. Plug the normal USB Connector (Large Connector) into the personal computer.

## 3. Plug the mini USB connector (Small Connector) into the USB connector.

“PC” will appear on the display when the USB cable is connected correctly.

### NOTE:


Follow the direction by instructor or service provider to use this function.

# CARE AND MAINTENANCE

## Clean the unit with a soft dry cloth.

- Do not use any abrasive or volatile cleaners.
- The unit is not waterproof. Do not wash it or touch it with wet hands.
- Never immerse the unit in water.
- Be careful so water does not get into the unit.

## Store the unit in a safe and dry location.

- Avoid the unit being exposed to direct sunlight, high temperatures, high humidity, water and/or dust.

## Do not subject the unit to strong shock, drop, or step on the unit.

## Remove the battery if the unit will not be used for three months or longer.

## Use the unit consistent with the instruction provided in this manual.

# TROUBLESHOOTING TIPS

<b>Problem</b>	<b>Probable cause</b>	<b>How to correct</b>
The Low Battery Indicator blinks.	The battery is low or worn.	Replace the battery with a new one (CR2032).
Nothing is displayed.	The polarities of battery (+ and -) are aligned in the wrong directions.	Insert the battery in correct alignment.
	The battery is worn.	Replace the battery with a new one (CR2032).
The displayed values are incorrect.	The unit is positioned incorrectly.	Refer to “Attaching the Pedometer” on page 16.
	Irregular movement. Inconsistent pace. Up and down movement.	
	The set values are wrong.	Check the weight and stride length settings.

# SPECIFICATIONS

Model:	<b>HJ-720ITCAN</b>
Power Supply:	3 VDC (1 lithium battery CR2032)
Measurement Range:	Steps: 0 to 99,999 steps Aerobic Steps: 0 to 99,999 steps Minutes for Aerobic Steps: 0 to 1,440 minutes Amount of Fat Burned: 0.0 to 199.9 g Calories: 0 to 99,999 kcal Distance: 0.00 to 999.99 km
Memory:	Previous 41 days (Most recent 7 days on display)
IT Function:	USB connector
Operating Temperature / Humidity:	-10°C to +40°C (14°F to 104°F) / 30% to 85% RH
Storage Temperature / Humidity / Air Pressure:	-20°C to +60°C (-4°F to 140°F) / 30% to 95% RH / 700 - 1600 hPa
Precision of Step Counting:	Within +/- 5% (by vibration testing machine)
Time Precision:	Within +/- 30 seconds of the average monthly deviation (under normal temperature)
Battery Life:	Approximately 6 months (when used for walking 10,000 steps a day)
External Dimensions:	Approximately 47 mm (w) × 73 mm (h) × 16 mm (d) (1 7/8" × 2 7/8" × 5/8")
Weight:	Approximately 37 g (battery included)
Contents:	Pedometer, Battery (CR2032), Screwdriver, Strap, Clip, Holder, and Instruction Manual
UPC:	0 73796 90314 5

- Specifications are subject to change without notice.

# FCC STATEMENT

## NOTE:

### POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the product and the receiver.
- Connect the product into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for Canada only)

This digital apparatus does not exceed the Class B limits for radio noise emissions from digital apparatus as set out in the interference-causing equipment standard entitled “Digital Apparatus”, ICES-003 of the Canadian Department of Communications.

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Changes or modifications not expressly approved by the party responsible for compliance could void the user’s authority to operate the equipment.

# WARRANTY INFORMATION

## Limited Warranty

Your Omron HJ-720ITCAN Pedometer, excluding the battery, is warranted to be free from defects in materials and workmanship appearing within 1 year from the date of purchase, when used in accordance with the instructions provided with the unit. The above warranties extend only to the original retail purchaser.

We will, at our option, repair or replace without charge any pedometer covered by the above warranties. Repair or replacement is our only responsibility and your only remedy under the above warranties.

To obtain warranty service contact Omron Healthcare's Customer Service by calling **1-800-634-4350** for the address of the repair location and the return shipping and handling fee. Information for warranty service is available on our website at [www.omronhealthcare.com](http://www.omronhealthcare.com).

Enclose the Proof of Purchase. Include a letter, with your name, address, phone number, and description of the specific problem. Pack the product carefully to prevent damage in transit. Because of possible loss in transit, we recommend insuring the product with return receipt requested.

**ALL IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE APPLICABLE WRITTEN WARRANTY ABOVE.** Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

**OMRON SHALL NOT BE LIABLE FOR LOSS OF USE OR ANY OTHER INCIDENTAL, CONSEQUENTIAL OR INDIRECT COSTS, EXPENSES OR DAMAGES.** Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above exclusions may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state.

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Copyright © 2006 OMRON HEALTHCARE, INC.